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## 📖 Enthält

17 mg x 2 Pflaster

MEDMETICS  
**REMESCAR**  
 INSTANT CORRECTIVE SKINCARE

**EYE BAGS &  
 DARK CIRCLES**  
 MICRONEEDLING  
 PATCHES

### 📄 Description

Remescar Eye Bags & Dark Circles Microneedling Patches are specially designed to reduce the appearance of eye bags and dark circles around the delicate eye area.

### Instructions for use

1. Cleanse and dry your skin. Peel off the patch from the film and remove the white protective sticker without touching the microstructure.
2. Apply the patch perpendicularly to the desired area and activate the patches by adding light pressure for few seconds.
3. Leave on for at least 4 hours or overnight. Carefully remove the patches.

Avoid using on the lower eyelid. Recommend to use twice a week. Clinically proven improvement of eye bags after 2 weeks.

### Precautions

- Keep out of reach of children.
- For external use only.
- Avoid contact with mucous membranes or the eyes. Rinse thoroughly with water in the event of contact with the eyes.
- Discontinue use if irritation or sensitization (redness, pain, etc.) occurs
- Only use on normal healthy, non-infected skin.
- Do not use if the packaging has been damaged.
- Not recommended for children under the age of 16.

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### Contenuto

17 mg x 2 cerotti

MEDMETICS  
**REMESCAR**  
INSTANT CORRECTIVE SKINCARE

**EYE BAGS &  
DARK CIRCLES**  
MICRONEEDLING  
PATCHES

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