remescar

OPPERVLAKKIGE
SPATADERS

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Ingrediënten:

Aqua, Caprylic/Capric Triglyceride (and) Titanium Dioxide (and) Polyhydroxystearic Acid (and) Aluminium Stearate (and) Alumina, Caprylic/Capric Triglyceride, Zinc Oxide (and) Caprylic/Capric Triglyceride (and) Polyhydroxystearic Acid (and) Isostearic Acid, Glycerin/Aqua/Palmaria palmata extract, Glyceryl Stearate/Behenyl Alcohol/Palmitic Acid/Stearic Acid/Lecithin/Lauryl Alcohol/Myristyl Alcohol/Cetyl Alcohol, Glycerin, Methylpropanediol (and) Caprylyl Glycol (and) Phenylpropanol, Squalane, Silica/Polyurethane-40/CI 61570, Tocopheryl Acetate, Aqua/Chlorella vulgaris extract/Phenoxyethanol, Sodium Ascorbyl Phosphate, PVP, Panthenol, Dimethicone, Aminomethyl propanol, Carbomer, Xanthan gum, Hydroxyethylcellulose, Sodium Hyaluronate, Hydrogenated Phosphatidylcholine, Citric Acid.

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Ingrédients:

Aqua, Caprylic/Capric Triglyceride (and) Titanium Dioxide (and) Polyhydroxystearic Acid (and) Aluminium Stearate (and) Alumina, Caprylic/Capric Triglyceride, Zinc Oxide (and) Caprylic/Capric Triglyceride (and) Polyhydroxystearic Acid (and) Isostearic Acid, Glycerin/Aqua/Palmaria palmata extract, Glycerol Stearate/Behenyl Alcohol/Palmitic Acid/Stearic Acid/Lecithin/Lauryl Alcohol/Myristyl Alcohol/Cetyl Alcohol, Glycerin, Methylpropanediol (and) Caprylyl Glycol (and) Phenylpropanol, Squalane, Silica/Polyurethane-40/CI 61570, Tocopheryl Acetate, Aqua/Chlorella vulgaris extract/Phenoxyethanol, Sodium Ascorbyl Phosphate, PVP, Panthenol, Dimethicone, Aminomethyl propanol, Carbomer, Xanthan gum, Hydroxyethylcellulose, Sodium Hyaluronate, Hydrogenated Phosphatidylcholine, Citric Acid.

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