

Description

Spider veins are small dilated blood vessels near the surface of the skin characterized as purple, blue or red lines. They are similar to varicose veins but are much thinner in diameter and not palpable.

Remescar Spider Veins Instant is an innovative product that has scientifically been formulated to decrease the appearance of superficial discolored, visual and structural skin imperfections due to spider veins. With only one application it will significantly decrease the appearance of spider veins by forming an instant film over the discolored skin area.

Try Remescar Spider Veins Instant and experience the self-assurance of a beautiful skin.

Intended use

Remescar Spider Veins Instant is designed to visually improve the skin affected by spider veins. It releases a transparent, thin and protective film which creates a barrier over the affected skin.

Instructions for use



Step 1: Wash your hands. The skin should be clean and dry.



Step 2: At first time use, open the tube and remove the protective cap.

Step 3: Squeeze out a small amount at a time to avoid wastage. Use more if you have a large area affected by spider veins.



Step 4: Massage the white cream gently in rotational motions on the affected skin. The cream will instantly conceal and decrease vascular imperfections. We recommend to use the product for 2-3 months for optimal results.



Step 5: Allow enough time for the cream to completely absorb into the skin and dry before dressing. Putting on clothes before the cream is fully absorbed can inhibit the cream's effectiveness and prevent it from creating a protective moisture barrier between your skin and the environment.



Step 6: Wash your hands to remove any remaining cream.

Discover Remescar's full range, including Remescar Crow's Feet, Remescar Eye Bags & Dark Circles, Remescar Stretch Marks and Remescar Scar Stick.

Precautions

- For external use only.
- Avoid contact with mucous membranes or the eyes.
- Keep out of the reach of children.
- In the case of vascular complications or varicose veins, consult a professional before use.
- Not to be used on dermatological conditions that disrupt the integrity of the skin (e.g. open wounds, third degree burns, acne and eczema).
- Use of the product is inadvisable if you have allergic reactions to one of the ingredients. If skin sensitization/irritation occurs, discontinue use and consult a professional.
- Leave at least 30 minutes after applying a moisturizing cream before you apply Remescar Spider Veins Instant on the same skin area, as otherwise the moisturizer could make Remescar Spider Veins Instant less effective.
- It's not advised to apply Remescar Spider Veins Instant over skin treatments.
- May stain clothes if not completely dry.
- Once opened, it's advised to use the product within 12 months.
- It's advisable to store the product below 30°C.
- In general, avoid extreme sun exposure of spider veins.

Ingredients

Aqua, Dibutyl Adipate, Dicaprylyl Carbonate, Isoamyl p-Methoxycinnamate, Glycerin, Methylpropanediol, Distarch Phosphate, Cellulose, Talc, Mannitol, Cetearyl Alcohol, Ethylhexyl Salicylate, Diethylamino Hydroxybenzoyl Hexyl Benzoate, Myristyl Myristate, Tocopherol Acetate, Glyceryl Stearate Citrate, Acrylates/Ethylhexyl Acrylate Copolymer, Hydrogenated Lecithin, Caprylyl Glycol, Polyglyceryl-3 Stearate, Cetearyl Glucoside, PEG-12 Dimethicone, Polymethyl Methacrylate, Tocopherol, Parfum, Sodium Hyaluronate, Xanthan Gum, Hydroxyethyl Acrylate/Sodium Acryloyldimethyl Taurate Copolymer, Phenylpropanol, Aluminium Hydroxide, Sodium Phytate, Chlorella Vulgaris Extract, Palmaria Palmata Extract, Laureth-21, Phenoxyethanol, Sorbitan Isostearate, Polysorbate 60, Glucose, Alcohol, CI 77491, CI 77492, CI 77499, CI 77891.

Manufacturer

Sylphar NV
Xavier De Cocklaan 42
9831 Deurle
Belgium
www.remescar.com

Distributor

Sylphar NV
Xavier De Cocklaan 42
9831 Deurle
Belgium
www.remescar.com

Contains

40 ml e