remescar

SPIDER VEINS

Description

Spider veins are small dilated blood vessels near the surface of the skin characterized as purple, blue or red lines. They are similar to varicose veins but are much thinner in diameter and not palpable.

Remescar Spider Veins is an innovative medical device designed for the prevention and treatment of superficial discolored, visual and structural imperfections (spider veins) by forming a long lasting breathable occlusive film over and around the affected skin. This thin and discrete barrier has several active properties.

The product will protect the skin from toxic aggressors and inhibit the vascular damage. Thanks to the phospholipid-induced spherical film, it will lead to an excellent integration of the cream into the skin lipid barrier.

In addition, Remescar Spider Veins will protect and provide structural support, so formation of new blood vessels will be prevented, the redness of the vascular lesions decreased and the microcirculation improved.

Thanks to these properties, Remescar Spider Veins will also reinforce and visually improve the skin surface.

Intended use

Remescar Spider Veins is designed for the prevention, management and treatment of skin affected by spider veins. It releases a transparent, thin and protective film which creates a barrier over the affected skin.

Instructions for use

The skin should be clean and dry. At first time use, open the tube and remove the protective cap. Wash your hands. Remescar Spider Veins should be applied gently on the affected skin. Massage the cream in until completely absorbed. Do not rub too hard as this will not increase the effectiveness. Apply Remescar Spider Veins twice a day for 2 to 3 months.

Precautions

- For external use only.
- Avoid contact with mucous membranes or the eyes.
- Keep out of the reach of children.
- In the case of vascular complications or varicose veins, consult a professional before use.
- Not to be used on dermatological conditions that disrupt the integrity of the skin (e.g. open

wounds, third degree burns, acne, eczema).

- Remescar Spider Veins is not sterile and must not be used as a wound dressing.
- Use of the product is inadvisable if you have allergic reactions to one of the ingredients. If skin sensitization/irritation occurs, discontinue use and consult a professional.
- Avoid using a topical cream and Remescar Spider Veins at the same time. It's advised to leave 30 minutes between both.
- It's not advised to apply Remescar Spider Veins over other skin treatments.
- May stain clothes if not completely dry.
- Do not use after the expiry date. Once opened, it's advised to use the product within 6 months.
- It's advisable to store the product below 30°C.
- In general, avoid extreme sun exposure of spider veins.

Ingredients

Aqua, Caprylic/Capric Triglyceride (and) Titanium Dioxide (and) Polyhydroxystearic Acid (and) Aluminium Stearate (and) Alumina, Caprylic/Capric Triglyceride, Zinc Oxide (and) Caprylic/Capric Triglyceride (and) Polyhydroxystearic Acid (and) Isostearic Acid, Glycerin/Aqua/Palmaria palmata extract, Glyceryl Stearate/Behenyl Alcohol/Palmitic Acid/Stearic Acid/Lecithin/Lauryl Alcohol/Myristyl Alcohol/Cetyl Alcohol, Glycerin, Methylpropanediol (and) Caprylyl Glycol (and) Phenylpropanol, Squalane, Silica/Polyurethane-40/CI 61570, Tocopheryl Acetate, Aqua/Chlorella vulgaris extract/Phenoxyethanol, Sodium Ascorbyl Phosphate, PVP, Panthenol, Dimethicone, Aminomethyl propanol, Carbomer, Xanthan gum, Hydroxyethylcellulose, Sodium Hyaluronate, Hydrogenated Phosphatidylcholine, Citric Acid.

Symbols

- = legal manufacturer
- = batch number
- = expiry date
- = read the instructions for use
- = keep away from sunlight
- = store below 30°C

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